



# NIDS® PEG IMMUNOGENICITY RAPID ASSAY

## Double Antigen Bridging Assay

Uses ANPEG-1 Monoclonal anti-PEG IgM Positive Control Ab

Optimized for Point of Care Patient Monitoring



- Lateral flow format
- Easy to use
- Handheld reader
- 40 minute turnaround time
- Limit of Detection: 125 ng/mL Ab
- No matrix effects observed
- Highly reproducible
- Extensively tested with human serum samples
- Detects Abs to PEGs of various lengths and with different linking groups

